


# 3D Nature Multilayer Zirconia Block



- 6 layers Multilayer color
- From 43% to 57% gradient for transparency
- Strength shows gradient from 700MPa to 1050MPa
- Suitable for anterior, crown and 7 unit bridges




Color Transmittance Strength


## Physical Characteristics

|                       |                            |
|-----------------------|----------------------------|
| Sintered density      | 6.07±0.03g/cm <sup>3</sup> |
| Bending Strength      | 700-1050MPa(Av.)           |
| Transmittance         | 57%-43%                    |
| Hardness              | 1200HV                     |
| Sintering temperature | Recommend 1500°C           |

## Color



## Adaptation systems



## Indications

