



*To Be More Natural
Explore,
what you really need*



4D PRO

Opening it, Full all your needs for denture indications!!!



Designed for full bridges



Tel : +86-13242096756
Email : info@yucera.com
Web : www.yucera.com

Shenzhen Yurucheng Dental Materials Co.,Ltd
No. 35, Zhuqing Road, Shijing Community, Shijing Street,
Pingshan District, Shenzhen, China

8 layers Multilayer color

15 layers shades gradient natural transparency



Layer 1	15%	57%	$\geq 700\text{MPa}$
Layer 2	10%	56%	$\geq 750\text{MPa}$
Layer 3	12%	54%	$\geq 820\text{MPa}$
Layer 4	12%	52%	$\geq 890\text{MPa}$
Layer 5	12%	50%	$\geq 960\text{MPa}$
Layer 6	12%	48%	$\geq 1030\text{MPa}$
Layer 7	12%	46%	$\geq 1100\text{MPa}$
Layer 8	15%	43%	$\geq 1200\text{MPa}$

Designed for full bridges

Color



A1 A2 A3 A3.5 A4



B1 B2 B3 B4



C1 C2 C3 C4



D2 D3 D4

BLEACH

OM1



OM2



OM3



Physical Characteristics

Sintered density	$6.07 \pm 0.03\text{g/cm}^3$
Bending Strength	$\geq 700\text{-}1200\text{MPa(Av)}$
Transmittance	43%-57%
Hardness	1200HV
Sintering temperature	Recommend 1500°C

Indications



Veneer



Crown



Anterior



Implant



14 unit
bridges